

# EM-900 WA SUPERVISED WIRELESS 40 PATIENT MONITOR STATION



**I recommend** the **EM-900WA Wanderer Alert** for patients, families and friends. The **NEW EM-900WA** Wanderer Alert monitors up to 40 patients at Long Range! It's portable, allowing continuous monitor while away from work station.

The **Wireless** pendant transmitters are worn by patients and continuously send signals to the Wanderer Alert. Whenever patients "Wander Away" or move outside areas or facilities, the signals are no longer received from the pendant transmitter. The Wanderer Alert detects the missing signals and **sounds** the alarm. The LCD displays the patient name, location, etc (the alpha-numeric information that you have entered) - received from the pendant transmitters. Our staff can quickly respond to the emergency.

A special feature **"Patient Locator"** is built into the Wanderer Alert and allows search and location in the event of missing or lost patients. Other advantages are "Call waiting" for multiple alarms and "Data Log" for analysis.


**Now... patients can be cared for without the need of confinement or physical restraints!**

We are **so happy with Wanderer Alert**...get it for your patients, you'll be glad you did!



Transmitter is worn on wrist or ankle with Locking Band

Pendant transmitter continuously sends the signal every 10 seconds. If signals are missing for 3 minutes, Wanderer Alert Sounds the Alarm!



EM-900WA is placed at Nurse station, desk or can be carried with you! Ni-Mh battery lasts up to 8 hours!




Wanderer Alert sounds the Alert and LCD displays information such as:

TX MISSING JOE AB

Event Name Location

You may enter any information up to 15 characters



"Wander Away." patients



"Missing" or "Lost" patients



"Walk Out" patients



Personnel or Items Leaving Facility

# Wireless 900 MHz Technology comes to the rescue of patients!

**Wanderer Alert Automatically Sounds the Alert and Notifies of who is missing,  
Even while care provider is roaming!**

EM-900 WA shown in  
re-charge cradle.  
Ni-Mh battery lasts up  
to 8 hours.



The EM-900 Wanderer Alert is **fast and easy** to set up and operate. Connect EM-900 cable to pendant transmitters. Enter information for each transmitter (ex: Event - Name - Location).

You can enter any information (up to 15 characters) that will be displayed on LCD screen. It will **Beep Alert and Display** entered information whenever a **transmitter is activated by patient going out of range. Unit is Portable!** Can be in hand, pocket or shoulder strap while roaming.



"Call Waiting" will retain up to 10 multiple incoming calls. Data Log retains calls for evaluation. Aux. out for Automatic Voice Dialers calls cellphones, pagers and telephones.

## How It Works:

1. The wireless pendant is attached to the patient with a non-removable band and **sends a continuous signal** to the Wanderer Alert unit every 10 seconds. The unit will Alarm if it doesn't hear a signal within 3 minutes.
2. When the patient attempts to wander away, the Wanderer Alert unit detects the missing signals and **sounds the Alert** at the EM-900. The LCD Displays the missing patient's name etc.
3. **"Patient Locator"** function will ask "Do you want to Search for the missing patient", if selected it will go into **"Search Mode"** allowing you to locate patient by detecting the signals from the missing patient's transmitter as you search.



**1 pendant per  
unit station only.  
(Max. 40)**

**The Wanderer Alert system is capable of monitoring  
both High Risk and Low Risk patients!**



### High Risk Patients \*patients that are high risk to wander away

For High Risk patients, the pendant is attached to ankle and will activate Wanderer Alert when patient wanders **up to 50 feet or less from unit.** Allowing patient to move throughout facility or living areas.

### Low Risk Patients \*patients that are low risk to wander.

With Low Risk patients, attach pendant to wrist. This allows them to move throughout living areas, facility, yard, etc. Wanderer Alert will activate when patient wanders **up to 500 feet or less from unit.**



## Order Information

Model	Description	
EM-900 WA	.. Complete with Wanderer Alert Unit, Charging Cradle, AC Adapter	
WPT-900	..... Wireless Pendant transmitter with extra Locking wrist-ankle bands	
LWA-900	..... Patient locator unit only with charging cradle, AC adapter	
AVD-45be	..... Automatic voice dialer	See Emergency Caller Product Index and Price list for all accessories for Products and Systems!

**Emergency Caller Products**



**Emergency Caller Products (T) 800-227-2474 (F) 800-762-7613**  
[www.emergencycaller.com](http://www.emergencycaller.com) [sales@emergencycaller.com](mailto:sales@emergencycaller.com)

**USP** United  
Security  
Products, Inc.

U.S.A. and Foreign Patents Issued and Pending  
Made In USA