
Patient Alert Fall Prevention System



Patient Alert is part of our economical Direct Wire Fall Prevention Program. Designed to operate with all Emergency Caller Sensor Pads it provides unmatched performance when protecting residents without the need of physical restraints.

- Selectable for resident attempts to get off or get on sensor pads.
- Sounds Alarm and Flashes a Strobe when activated to remind resident to return to safe location or position, and notify staff or caretaker.
- Instant or Delayed Activation Models.
- Selectable alarm 5 seconds or continuous until turned off by caregiver.
- Low battery alert uses long life 9 volt alkaline battery.
- Nurse Station Output Jack

ORDERING INFORMATION:

PA

Caller Products


Ready to Go systems or customize your own!

- PA-1 Patient Alert Unit, 1" x 30" strip sensor pad, accessories
- PA-2 Patient Alert Unit, 9" x 16" Wheelchair/Seat sensor pad, 4 ft. cable
- PA-3 Patient Alert Unit, 14" x 24" Bed/Floor sensor pad, accessories
- PA-4 Patient Alert Unit, 18" x 24" Bed/Floor sensor pad, accessories
- PA-5 Patient Alert Unit, 24" x 30" Bed/Floor sensor pad, accessories
- NCC-8 1/4" phone jack, 8 ft cable adapts Patient Alert to Nurse Call Station
- XC-25 Extension cable allows Patient Alert additional 25 ft distance from sensor pad
- NCA-22 "Y" Adaptor -- adapts PA-(X) units and Nurse Call Button to Nurse Call Station: 1/4 mono
 - Additional Patient Alert Unit Only 2 second delay
- PA-ND Additional Patient Alert Unit Only No Delay (instant activation) *above systems may also be ordered with PA-ND no delay unit

* Soft Touch cover included on PA-2, PA-3 and PA-4 systems, as well as P-6-22, P-3-22 & P-4-22 pads sold separately.

Please see Sensor Pad page for additional Sizes

Emergency Caller Products (I) 800-227-2474 (F) 800-762-7613





Emergency

PATIENT ALE